

ELK Biotechnology

mStrawberry-Tag Mouse mAb

Catalog NO.: EM1164 For research use only.

Overview

Product name mStrawberry-Tag Mouse Monoclonal antibody

Source Mouse

Applications WB

Species reactivity Human Rat Mouse

Recommended dilutions WesternBlot:1/2000-5000

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species N/A

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

lsotype lgG1

Clonality Monoclonal

Concentration 1 mg/ml

Observed band N/A

GenelD (Human) N/A

Human Swiss-Prot No. N/A

Cellular localization N/A

Alternative Names mStrawberry antibody

Background mStrawberry is a bright monomeric red fluorescent protein which was

developed in Dr. Roger Tsien's lab by directed mutagenesis of mRFP.



Western blot analysis of mStrawberry recombinant protein with mStrawberry (Mix-mA $^{\text{TM}}$) Mouse mAb EM1164 diluted at:5000.