



# ELK Biotechnology

ECFP-Tag Mouse mAb

Catalog NO.: EM1178

For research use only.

## Overview

---

Product name	ECFP-Tag Mouse Monoclonal antibody
Source	Mouse
Applications	<b>WB</b>
Species reactivity	<b>N/A</b>
Recommended dilutions	<b>WesternBlot:1/2000-5000</b> <b>NOTE: Optimal dilutions should be determined by the end user.</b>
Immunogen	Recombinant Protein
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1 mg/ml
Observed band	<b>N/A</b>
GeneID (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	Cyan fluorescent protein Enhanced Cyan fluorescent protein
Background	Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish <i>Aequorea victoria</i> .

1

2

Western blot analysis of Recombinant ECFP Protein with ECFP Mouse mAb EM1178 diluted at)1:5000 2)1:10000

