

Recombinant Human PTH

Catalog #	EPT146
Expression Host	E.coli
DESCRIPTION	Recombinant Human Parathyroid Hormone is
	produced by our E.coli expression system and the
	target gene encoding Ser32-Gln115 is expressed.
Accession	P01270
Synonyms	Parathyroid Hormone; PTH; Parathormone; Parathyrin
Mol Mass	9 KDa
AP Mol Mass	14 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL
	test.
FORMULATION	Lyophilized from a 0.2 μ m filtered solution of 10mM
	HAc-NaAc, 150mM NaCl, 5% Mannitol, pH 4.0.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by
	vortex or pipetting.

It is not recommended to reconstitute to a



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com 23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q



concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water.

Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING The product is shipped at ambient temperature.Upon receipt, store it immediately at the temperature listed below.

STORAGELyophilized protein should be stored at < -20 ° C,
though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at 4-7°C
for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

BACKGROUND Parathyroid hormone is the most important endocrine regulator of calcium and phosphorus concentration in extracellular fluid. This hormone is secreted from cells of the parathyroid glands and finds its major target cells in bone and kidney. Another hormone, parathyroid hormone-related protein, binds to the same receptor as parathyroid hormone and has major effects on development. Like most other protein hormones, parathyroid hormone is synthesized as a



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



preprohormone. After intracellular processing, the mature hormone is packaged within the Golgi into secretory vesicles, the secreted into blood by exocytosis. Parathyroid hormone is secreted as a linear protein of 84 amino acids.



SDS-PAGE



+86-27-59760950 ELKbio@ELKbiotech.com

ch.com www.e

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C