

## Recombinant Mouse B7-H4 (C-Fc)

Catalog #	EPT204
Expression Host	Human Cells
DESCRIPTION	Recombinant Mouse B7 Homolog 4 is produced by
	our Mammalian expression system and the target
	gene encoding Phe29-Pro258 is expressed with a Fc
	tag at the C-terminus.
Accession	AAH32925.1
Synonyms	V-set domain containing T-cell activation inhibitor 1;
	B7 homolog 4; Protein B7S1; B7-H4; VTCN1
Mol Mass	52.2 KDa
AP Mol Mass	70-90 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ $\mu$ g (1 EU/ $\mu$ g) as determined by LAL
	test.
FORMULATION	Lyophilized from a 0.2 $\mu m$ filtered solution of PBS, pH
	7.4.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com 23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



vortex or pipetting.

It is not recommended to reconstitute to a concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING The product is shipped at ambient temperature.Upon receipt, store it immediately at the temperature listed below.

STORAGELyophilized protein should be stored at < -20 ° C,<br/>though stable at room temperature for 3 weeks.<br/>Reconstituted protein solution can be stored at 4-7°C<br/>for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

BACKGROUND Mouse V-set domain-containing T-cell activation inhibitor 1/VTCN1/B7-H4 is glycosylated member of the B7 family of immune co - stimulatory proteins. B7-H4 consists of extracellular domain (ECD) with one Ig-like V-set domain and one Ig-like C2-set domain. It is widely expressed, including in kidney, liver, lung, pancreas, placenta, prostate, spleen, testis and



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



thymus. B7-H4 negatively regulates T-cell-mediated immune response by inhibiting T-cell activation, proliferation, cytokine production and development of cytotoxicity. When expressed on the cell surface of tumor macrophages, plays an important role, together with regulatory T-cells (Treg), in the suppression of tumor-associated antigen-specific T-cell immunity. It also involved promoting epithelial cell in transformation.



## **SDS-PAGE**



ELKbio@ELKbiotech.com +86-27-59760950

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q