

Recombinant Human VEGF-D (C-6His)

Catalog #	EPT213
Expression Host	Human Cells
DESCRIPTION	Recombinant Human Vascular Endothelial Growth
	Factor D is produced by our Mammalian expression
	system and the target gene encoding Phe93-Ser201 is
	expressed with a 6His tag at the C-terminus.
Accession	O43915
Synonyms	Vascular Endothelial Growth Factor D; VEGF-D;
	c-Fos-Induced Growth Factor; FIGF; VEGFD
Mol Mass	13 KDa
AP Mol Mass	18 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL
	test.
FORMULATION	Lyophilized from a 0.2 μm filtered solution of PBS, pH
	7.4.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



vortex or pipetting.

It is not recommended to reconstitute to a concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING The product is shipped at ambient temperature.Upon receipt, store it immediately at the temperature listed below.

STORAGELyophilized protein should be stored at < -20 ° C,
though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at 4-7 °C
for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

www.elkbiotech.com

BACKGROUND Vascular endothelial growth factor D (VEGF-D) is a member of the platelet-derived growth factor/vascular endothelial growth factor (PDGF/VEGF) family. It is highly expressed in lung, heart, small intestine and fetal lung, and at lower levels in skeletal muscle, colon, and pancreas. VEGF-D is growth factor active in angiogenesis, lymphangiogenesis and



+86-27-59760950

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q

ELKbio@ELKbiotech.com



endothelial cell growth, stimulating their proliferation and migration and also has effects on the permeability of blood vessels. It may function in the formation of the venous and lymphatic vascular systems during embryogenesis, and also in the maintenance of differentiated lymphatic endothelium in adults. It undergoes complex proteolytic maturation, а generating multiple processed forms that bind and activate VEGFR-2 and VEGFR-3 receptors.



SDS-PAGE



ELKbio@ELKbiotech.com +86-27-59760950

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q