

Recombinant Mouse GH

Catalog #	EPT217
Expression Host	E.coli
DESCRIPTION	Recombinant Mouse Growth Hormone is produced by
	our E.coli expression system and the target gene
	encoding Phe27-Phe216 is expressed.
Accession	P06880
Synonyms	Somatotropin; Growth Hormone; Gh1; Gh
Mol Mass	22 KDa
AP Mol Mass	20 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL
	test.
FORMULATION	Lyophilized from a 0.2 μm filtered solution of 20mM
	Tris-HCl, 150mM NaCl, pH 8.0.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by
	vortex or pipetting.

It is not recommended to reconstitute to a



ELKbio@ELKbiotech.com www.elkbiotech.com +86-27-59760950 23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize

freeze-thaw cycles.

SHIPPING The product is shipped at ambient temperature.Upon receipt, store it immediately at the temperature listed below.

STORAGELyophilized protein should be stored at < -20 ° C,
though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at 4-7°C
for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

BACKGROUND Somatotropin(GH) is a member of the somatotropin/prolactin family of hormones which play an important role in growth control. Its major role in stimulating body growth is to stimulate the liver and other tissues to secrete IGF-1. GH stimulates both the differentiation and proliferation of myoblasts. It also stimulates amino acid uptake and protein synthesis in muscle and other tissues.



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com 23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C







+86-27-59760950 ELKbio@E

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C