

## **Biotinylated Human Frizzled-8**

## (C-6His-Avi)

Catalog #	EPT227
Expression Host	Human Cells
DESCRIPTION	Biotinylated Recombinant Human Frizzled-8 is
	produced by our Mammalian expression system and
	the target gene encoding Ala28-Pro172 is expressed
	with a 6His, Avi tag at the C-terminus.
Accession	Q9H461
Synonyms	frizzled 8; frizzled family receptor 8; Frizzled8;
	Frizzled-8; FZ-8; FZD8; hFZ8
Mol Mass	19.1 KDa
AP Mol Mass	25-35 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ $\mu$ g (1 EU/ $\mu$ g) as determined by LAL
	test.
FORMULATION	Lyophilized from a 0.2 $\mu\text{m}$ filtered solution of PBS, pH
	7.4.



www.elkbiotech.com +86-27-59760950 23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q

ELKbio@ELKbiotech.com



RECONSTITUTION Always centrifuge tubes before opening.Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. SHIPPING The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below. Lyophilized protein should be stored at < -20 °C, STORAGE though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at  $< -20^{\circ}$ C for 3 months. BACKGROUND Frizzled-8 is one of at least ten seven-transmembrane (7TM) glycoproteins of the Frizzled family of Wnt receptors. Frizzled proteins are thought to be G-protein-coupled. Wnt engagement, with low density lipoprotein receptor-related proteins LRP-5 or LRP-6 acting as co-receptors, stabilizes beta -catenin



+86-27-59760950

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q

www.elkbiotech.com

ELKbio@ELKbiotech.com



and promotes gene transcription that is important in development and tissue maintenance. Component of the Wnt-Fzd-LRP5-LRP6 complex that triggers beta-catenin signaling through inducing aggregation of receptor-ligand complexes into ribosome-sized signalosomes. The beta-catenin canonical signaling pathway leads to the activation of disheveled proteins, inhibition of GSK-3 kinase, nuclear accumulation of beta-catenin and activation of Wnt target genes. These ligands bind the extracellular CRD of Frizzled-8, blocking Wnt binding. The recombinant Frizzled-8 CRD has also been used to block Wnt signaling and inhibit growth of teratocarcinomas.



## SDS-PAGE



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C