

## Recombinant Human VDB (C-6His)

Catalog #	EPT229
Expression Host	Human Cells
DESCRIPTION	Recombinant Human Vitamin D-Binding Protein is
	produced by our Mammalian expression system and
	the target gene encoding Leu17-Leu474 is expressed
	with a 6His tag at the C-terminus.
Accession	P02774
Synonyms	Vitamin D-Binding Protein; DBP; VDB; Gc-Globulin;
	Group-Specific Component; GC
Mol Mass	52.3 KDa
AP Mol Mass	53 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ $\mu$ g (1 EU/ $\mu$ g) as determined by LAL
	test.
FORMULATION	Lyophilized from a 0.2 $\mu\text{m}$ filtered solution of PBS, pH
	7.4.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



vortex or pipetting.

It is not recommended to reconstitute to a concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING The product is shipped at ambient temperature.Upon receipt, store it immediately at the temperature listed below.

STORAGELyophilized protein should be stored at < -20 ° C,<br/>though stable at room temperature for 3 weeks.<br/>Reconstituted protein solution can be stored at 4-7°C<br/>for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

**BACKGROUND** Vitamin D-Binding Protein (DBP) is a member of the ALB/AFP/VDB family. DBP is a secreted protein and contains three albumin domains. The primary structure contains 28 cysteine residues forming multiple disulfide bonds. DBP acts as a multifunctional protein found in plasma, ascitic fluid, cerebrospinal fluid, and urine and on the surface of many cell types.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



DBP binds to vitamin D and its plasma metabolites and transports them to target tissues. DBP associates with membrane-bound immunoglobulin on the surface of B-lymphocytes and with IgG Fc receptor on the membranes of T-lymphocytes.



## **SDS-PAGE**



+86-27-59760950 ELKbio@ELKb

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C