

## Recombinant Human IL-2 (C-6His)

Catalog #	EPT246
Expression Host	Human Cells
DESCRIPTION	Recombinant Human Interleukin-2 is produced by our
	Mammalian expression system and the target gene
	encoding Ala21-Thr153 is expressed with a 6His tag at
	the C-terminus.
Accession	P60568
Synonyms	Interleukin-2; IL-2; T-Cell Growth Factor; TCGF;
	Aldesleukin; IL2
Mol Mass	16.4 KDa
AP Mol Mass	14-18 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.001 ng/ $\mu$ g (0.01 EU/ $\mu$ g) as determined by
	LAL test.
FORMULATION	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM
	PB, 150mM NaCl, pH 7.4.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com 23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q



vortex or pipetting.

It is not recommended to reconstitute to a concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING The product is shipped at ambient temperature.Upon receipt, store it immediately at the temperature listed below.

STORAGELyophilized protein should be stored at < -20 ° C,<br/>though stable at room temperature for 3 weeks.<br/>Reconstituted protein solution can be stored at 4-7°C<br/>for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

www.elkbiotech.com

**BACKGROUND** Interleukin-2(IL-2) is an interleukin, a type of cytokine signaling molecule in the immune system, belongs to the IL-2 family. It is a powerful immunoregulatory lymphokine produced by T-cells in response to antigenic or mitogenic stimulation. IL-2/IL-2R signaling is required for T-cell proliferation and other fundamental functions that are essential for the



+86-27-59760950

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q

ELKbio@ELKbiotech.com



immune response. IL-2 stimulates growth and differentiation of B-cells, NK cells, lymphokine-activated killer cells, monocytes, macrophages and oligodendrocytes.



## **SDS-PAGE**



+86-27-59760950 ELKbio@ELKbiotech.com

otech.com www

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C