

## Recombinant Human ST2 (C-6His)

Catalog #	EPT263
Expression Host	Human Cells
DESCRIPTION	Recombinant Human Interleukin-1 Receptor-like 1 is
	produced by our Mammalian expression system and
	the target gene encoding Lys19-Phe328 is expressed
	with a 6His tag at the C-terminus.
Accession	Q01638-2
Synonyms	Interleukin-1 receptor-like 1;Protein ST2;DER4; ST2; T1
Mol Mass	36 KDa
AP Mol Mass	50-85 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ $\mu$ g (1 EU/ $\mu$ g) as determined by LAL
	test.
FORMULATION	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM
	PB, 150mM NaCl, pH 7.4.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by
	vortex or pipetting.



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com 23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



It is not recommended to reconstitute to a concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPINGThe product is shipped at ambient temperature.Upon receipt, store it immediately at the temperaturelisted below.

STORAGELyophilized protein should be stored at < -20 ° C,<br/>though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7°C<br/>for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

BACKGROUND Interleukin 1 receptor-like 1(IL1RL1) is a member of the interleukin-1 receptor family, Contains 3 Ig-like C2-type domains and 1 TIR domain. It is highly expressed in kidney, lung, placenta, stomach, skeletal muscle, colon and small intestine. IL1RL1 is a receptor for interleukin-33, its stimulation recruits MYD88, IRAK1, IRAK4, and TRAF6, followed by phosphorylation MAPK3/ERK1 and/or of



+86-27-59760950

ELKbio@ELKbiotech.com

iotech.com www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.G



MAPK1/ERK2, MAPK14, and MAPK8. IL1RL1 may possibly be involved in helper T-cell function. Soluble IL1RL1 also acts as a negative regulator of Th2 cytokine production, it directly implicated in the progression of cardiac disease.



## **SDS-PAGE**



+86-27-59760950 ELKbio@ELKbiotech.com

h.com www.elk

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C