

## Recombinant Human MetAP1D (N,

## C-6His)

Catalog #	EPT270
Expression Host	E.coli
DESCRIPTION	Recombinant Human Methionine Aminopeptidase 1D
	Mitochondrial is produced by our E.coli expression
	system and the target gene encoding Arg44-Ala335 is
	expressed with a 6His tag at the N-terminus, 6His tag
	at the C-terminus.
Accession	Q6UB28
Synonyms	Methionine Aminopeptidase 1D Mitochondrial;
	Methionyl Aminopeptidase Type 1D Mitochondrial;
	METAP1D; MAP1D
Mol Mass	35.4 KDa
AP Mol Mass	35-40 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ $\mu$ g (1 EU/ $\mu$ g) as determined by LAL
	test.



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



FORMULATION	Supplied as a 0.2 $\mu$ m filtered solution of 20mM
Market Bark	Tris-HCl, pH 8.0.
RECONSTITUTION	
SHIPPING	The product is shipped on dry ice/polar packs.
	Upon receipt, store it immediately at the temperature
	listed below.
STORAGE	Store at $\leq$ -70°C, stable for 6 months after receipt.
	Store at $\leq$ -70 °C, stable for 3 months under sterile
	conditions after opening.
	Please minimize freeze-thaw cycles.
BACKGROUND	Methionine Aminopeptidase 1D (METAP1D) is a
	mitochondrion protein that belongs to the peptidase
	M24A family. METAP1D is overexpressed at the
	protein level in colon cancer cell lines and colon
	tumors as compared to normal tissues. N-terminal
	methionine removal is an important cellular process
	required for proper biological activity, subcellular
	localization, and eventual degradation of many
	proteins. METAP1D is also active with zinc, manganese
	or divalent ions. It may also play an important role in
	colon tumorigenesis.



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com 23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





## SDS-PAGE



+86-27-59760950 ELKbio@El

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C