

Recombinant Mouse C3a

Catalog #	EPT296
Expression Host	E.coli
DESCRIPTION	Recombinant Mouse Complement Component C3a
	is produced by our E.coli expression system and the
	target gene encoding Ser671-Arg748 is expressed.
Accession	P01027
Synonyms	Complement Component C3a; Anaphylatoxin; C3a
Mol Mass	9.2 KDa
AP Mol Mass	9 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL
	test.
FORMULATION	Lyophilized from a 0.2 μm filtered solution of PBS, pH
	7.4.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by
	vortex or pipetting.
	It is not accounted to the state of



It is not recommended to reconstitute to a



concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING The product is shipped at ambient temperature.Upon receipt, store it immediately at the temperature listed below.

STORAGELyophilized protein should be stored at < -20 ° C,
though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at 4-7°C
for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

BACKGROUND Complement is defined as key part of innate immunity and as the first line of defense in the fight against invading pathogens. Complement 3 (C3) is the most abundant component of the complement cascade and the convergent point for all three major complement activation pathways: namely classical, alternative and mannose-binding lectin pathways. Complement activation leads to the formation of the C3 convertase, which cleaves C3 into the key effector molecules, C3a



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkk

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



(anaphylatoxin) and C3b (opsonin) which then drive microbe removal. By binding to C3a receptor (C3aR), C3a exhibits potent anaphylatoxin activity, including increased vascular permeability, triggering degranulation of mast cells, inflammation, and activating leukocytes.



SDS-PAGE



+86-27-59760950 ELKbio@ELKbiotech.com

piotech.com wv

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C